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PORTABLE SURGICAL ASPIRATOR USER MANUAL

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Please read this manual carefully before use

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SERİ NO

CONTENTS

SUBJECTS	PAGE NO
1. INTENDED USE	5
2. CAUTIONS	5
3. DEVICE DESCRIPTION	7
4. ACCESSORIES	11
5. TECHNICAL FEATURES	11
6. SETUP AND CONTROLS	13
7. APPLICATION	15
7.1. Setting Vacuum Levels	15
7.2. Evacuation of the Jars	15
8. SAFETY SYSTEMS	17
9. REQUIREMENTS AT USE	17
10. MAINTANENCE AND CLEANING	19
11. TROUBLESHOOTING	19
12. HANDLING AND STORAGE	21
13. WARRANTY TERMS	21
14. COPYRIGHT	23
15. CERTIFICATES	23
16. MANUFACTURER COMPANY INFORMATION	23

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1. INTENDED USE


LIFETIME Surgical Aspirator Devices are used for picking up the fluids and particles in surgical stage to a separate collection jar by the negative pressure generated with the vacuum pump.

LIFETIME Surgical Aspirator Devices are used in; polyclinics, clinics and by patients at home.

2. CAUTIONS

- Please make sure that you have read and understand this user manual along with its cautions and descriptions before setup.
- The person buying this device is responsible for the content of this manual and informing instructions to the relevant personnel.
- Do not use the device for any other purpose.
- The device should be operated by trained personnel.
- High Vacuum and High Flow generating device can be used in continuous or discontinuous modes.
- Perform periodic maintenance of the device on regular basis.
- Replace the filter after septic cases.
- If hydrophobic filter is used in the device, fluid contact will stop the aspiration system. In this case, please change the filter.
- Precautions should be taken at the power input against sudden voltage drops. Voltage regulator is recommended.

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
3.DEVICE DESCRIPTION


ON / OFF BUTTON : It is used to start up the vacuum pump.

VACUUM PORT : It is used for vacuum.


VACUUM CONTROL KNOB : It is used to set the vacuum level.

VACUUM GAGE: It is used to show the vacuum flow when the device is working.

VACUUM JAR : It is used to collect the aspirated fluids.

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


BAĞLANTI
APARATI

KATETER


Please make the bacterial filter or hydrophobic filter connection as shown.

Hose connections on the jar ;
1.VACUUM : Connect the hose with the bacterial filter or hydrophobic filter.
2. PATIENT : Connect the hose which the patient will use.

Aspiration Hose –
Connection Apparatus–
Catheter connection as
shown in the figure

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6.SETUP AND CONTROLS

Open the box and take out the device and other accessories from the box.

Check for any damages that can be occurred during transport. Please refer to the related sales point.

Connect one end of the hose with a bacterial filter or hydrophobic filter to the vacuum port and the other end to the port marked as 'vacuum' on the jar lid. Connect the long patient hose to the 'patient' port of the cover. Connect the aspiration hose, connection apparatus and catheter.

Please check the following controls in the first setup and before each use;

- Check to see if there is any cracks on the jar lid that will lead to an air leakage.
- Check the jar lid and overflow protection system.
- Make sure the jar lid seal is connected to the jar lid and jar lid is correctly connected to the jar.
- Make sure the jar is correctly mounted to the device.
- Check if there is any wearing, opening etc. on the aspiration hose.
- Check for any damages on the electrical cable. Do not use damaged electrical cables.
- Turn on the ON/OFF and make sure the right or left vacuum buttons illuminate.
- Press Start button to run the vacuum pump.
- Check the vacuum gauge by blocking the suction port.

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7. APPLICATION

Plug-in the electric cord of the device.

Turn on the OFF/ON button.

After switching ON / OFF button, the device starts running.

To stop the aspiration, press ON/OFF button again.

7.1. Setting the vacuum levels.

The vacuum adjustment system is available in the device. This helps you to use the device in the desired vacuum value. To make the adjustment please run the device and the vacuum pump will be activated. Please read vacuum gauge after closing the tip of the suction port. Adjust the setting to desired value using the adjusting knob.

7.2. Evacuation of Collection Jars

During the use of the device please control the waste accumulated in the collection jar in regularly. The collected waste level should not exceed the "Maximum" level mark on the jar. To evacuate the collection jar;

Turn off the device and remove the input-output hoses from the collection jar lid. Take off the collection jar from the plastic holder. Open the lid of the collection jar after taking it outside the clinic. Consider the collected waste liquid inside the collection jar as "contaminated medical waste" and always evacuate carefully outside clinical zone. Please subject the discharged collection jar disinfection treatment. Cover the lid and be sure that it has been closed tightly otherwise, leakage will occur between the collection jar and it will reduce device's vacuum performance.

This instruction should strictly be followed to avoid overflow of the waste liquid and elope into the pump during evacuation process. Otherwise, as the performance of the device will reduce but this will cause pump breakdown which will shorten the life of the device.

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8. SAFETY SYSTEMS

1. Float System: It's located under the collection jar lid. If the accumulated waste exceeds the maximum level in the jar and the evacuation does not take place then the float system actuated and protect the vacuum pump from waste. When this system is activated, the device stops working and the jar do not collect excess waste.

2. Bacteria Filter: Mounted on the hose between the device and the collection jar. Filter the air vacuumed from collection jar towards the device and helps to dispose clean air to the environment. At the end of its life the filter is blocked and the devices can't perform suction. Filter should be replaced with a new one.

3. Hydrophobic Filter: This is an optional system mounted on the hose in between the device and the collection jar. It filters the air vacuumed from collection jar toward the device and helps to dispose clear air to the environment. At the end of its life the filter is blocked and the devices can't perform suction. Filter should be replaced with a new one. It also prevents escape of waste to the device from the collection jar.

9. REQUISITES DURING APPLICATION

To avoid misuse and applications which may harm the patient, user or device guidelines should be strictly followed.

Do not damage the power cord. Please control robustness of the power cord before use each time. Do not use damaged power cords.

Do not use too long (more than 2 meters) aspiration tube. Increase of aspiration tube length may cause to low vacuum level.

Avoid applying any weight or force on the aspiration hose that will cause the hose to become clogged.

Please use clean aspiration hose and collection jar each time. In case of emergency please keep enough collection jars and aspiration hose in reserve. At the end of the use of the device the OFF/ON switch located at the back of the device should be turned off.

Do not use the device near explosives or in Oxygen-rich environment. Otherwise, it may elicit the danger of explosion or fire.

Perform the periodic maintenance of the device regularly by authorized personnel.

Do not run the device once the waste level reaches to the maximum before discharge the collection jar.

Do not use the device for the purpose other than described in this manual.

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10. MAINTENANCE AND CLEANING

Please be sure that the device has been switched off and the power cord has been unplugged before proceed for maintenance and cleaning of the device. Wear protective gloves during this procedure.

10.1 CLEANING OF THE DEVICE

Each time before using the device disinfect the outer surface of the device, collection jars and lids using medical disinfectant. The disinfectant should be used as per its own operating instructions.

10.2. STERILIZATION

The collection jars of the device can be autoclaved at 120°C after discharging the collected waste and cleaning the jars.

10.3. CHANGE OF FILTER

Replace the bacteria filters and the hydrophobic filters used in the device with new ones when soiled or obstructing the vacuum process.

11. TROUBLESHOOTING

FAULT	CAUSE	SOLUTION
Motor is not running	No power in the motor Power cord is disconnected Faulty motor	Check the fuse of the device check the power cord and the plug contact authorized technic service
Motor running, low vacuum value	Dirty Filters Filters has been clogged Collection jar filled up Lids of the jars not placed properly Faulty vacuum meter	Replace the filters Discharge the waste from the jar Check the lids Contact authorized technic service
Abnormal sound from the device	Defective motor Worn-out motor vibration dampers Waste entered into the motor	Contact authorized technic service

For SA03PT model, unplug the power cord of the device from the plug socket and pull off the fuse cover located below the input power socket using a screwdriver. Remove the defective fuse and replace with a new one. Cover the fuse by pushing it off.

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12. HANDLING AND STORAGE

The device and accessories must be transported in its original box. If you have no original box please make a suitable package using Styrofoam and carton box to prevent damage during transport of the device. The device should not be turned upside down during transport. Please act in accordance with the instructions contained on the box of the device during transportation. Store the device at normal room condition.

13. WARRANTY TERMS

1. Warranty period(2years) starts from the date of delivery of the goods to customers.
2. If the device fails within warranty period, the time spent in repairing is added to the warranty period. The repair period is maximum 30 working days
3. If you have a complaint, please send the complete device to your specialist dealer or send it to us in the original packaging, postage paid including warranty certificate inside with the specialist dealers stamp.
4. Damages arising from improper use of the device are not covered by this warranty
5. This warranty will be invalidated if the device is interfered with by unauthorised persons
6. Malfunction from not following the instructions in the manual are not covered by the warranty.
7. This warranty is valid only if all the conditions stated herein are read and accepted.
8. Consumable parts are not covered under warranty

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ASO 1.OSB BABÜRŞAH CD NO : 17 SİNCAN ANKARA

14. COPYRIGHT

This document has been issued by ELMASLAR TIBBI CIHAZ IMALAT A.S. and reserves the right to make changes to product specifications. No part of this publication may be reproduced or copied without written permission of the company.

15. CERTIFICATES

ISO 9001 , ISO 13485 , CE

16. MANUFACTURER COMPANY INFORMATION

ELMASLAR A.Ş.

ASO 1.OSB BABÜRSAH CD NO : 17 SINCAN ANKARA / TURKEY

ELMASLAR A.S.

ASO 1.OSB BABURSAH CD. NO: 17 SINCAN / ANKARA / TURKEY

Tel: +90 312 394 80 01 Fax: +90 312 394 80 04

info@elmaslarmedikal.com.tr • elmaslarmedical@gmail.com

www.elmaslarmedikal.com.tr