MODEL: Q5

**OKSİJEN KONSANTRATÖRÜ KULLANIM KILAVUZU
OXYGEN CONCENTRATOR USER GUIDE**

| CE 1984

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Please read this manual thoroughly before use of the Device

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SUBJECTS

PAGE NO

1. GENERAL INFORMATION	7
2. SAFETY GUIDELINES	7
3. INSTALLATION / USE	11
4. PREVENTIVE MAINTENANCE	13
5. TROUBLESHOOTING	17
6. AIRFLOW DIAGRAMS	19
7. TECHNICAL CHARACTERISTICS	21
8. LIST OF ACCESSORIES	21
9. ELECTROMAGNETIC IMMUNITY	23
10. GUARANTEE (TERMS & CONDITIONS)	27
11. CERTIFICATES	27
12. MANUFACTURER - AUTHORIZED TECHNICAL SERVICES	27

GLOSSARY OF SYMBOLS

	CAUTION CONSULT ACCOMPANYING DOCUMENTS		CLASS II EQUIPMENT		CONFORMITE EUROPE		NOT MOVED WITH HOOK
	DATE OF MANUFACTURE		PROTECTION AGAINST VERTICALLY FALLING DROPS		THIS SIDE UP		SHOULD BE DISPOSED ACCORDING TO WEEE DIRECTIVES
	MANUFACTURE		TYPE B		FRAGILE, HANDLE WITH CARE		KEEP AWAY FROM FIRE
	SERIAL NUMBER		TYPE BF		KEEP DRY		NO SMOKING

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1. GENERAL INFORMATION

Before operating the device, make sure that you have carefully read and understood all the topics written in this manual.

Intended Use of the Device: It is used to provide oxygen therapy to the patients that oxygen therapy is required (such as heart, lung and other respiratory diseases) whom determined and recommended by the doctor.

Working Principle: This oxygen concentrator is designed for personal use at home by providing the separated oxygen from the air in the room. This device obtains 90% pure oxygen from the environment at the rate of 21% oxygen by using PSA technology. Thanks to the cannula attached to the nose, it allows the patient to receive highly concentrated oxygen directly.

Side Effects: As with any treatment, oxygen therapy has risks and unwanted side effects. The most important risks in oxygen treatment are accidents and fires related to the transportation, filling and use of oxygen devices. Since oxygen is a flammable and explosive substance, the oxygen source should be placed away from heat generating systems and the patient should never smoke during oxygen therapy. Potential side effects of exposure to high concentrations of oxygen are changes in normal physiological functions, oxygen-related tissue damage / oxygen toxicity, and carbon dioxide accumulation. Carbon dioxide accumulation is mostly seen as a result of using high levels of oxygen, and it is usually manifested by drowsiness in patients.

Patients Need Oxygen Therapy: Patients with Chronic Obstructive Pulmonary, Pulmonary Hypertension, Congestive Heart Failure, Cystic fibrosis, Lung cancer, Alpha 1 and similar patients.

Patient Population: As of 2015, COPD affects 174.5 million people in the world and this number corresponds to 2.4% of the world population.

It occurs especially in people over the age of 40. It is seen at almost the same rates in men and women. COPD, which killed 2.4 million people in 1990, was the cause of death for 3.2 million people in 2015.

General Features

It has a Safety and Solid body

LCD display showing the total usage time and fault notifications

Electrical cut-off warning for safety against high pressure, low pressure & high temperature.

Thermo protected compressor,

Light and easy to carry.

Contraindications: Untreated pneumothorax patients, patients using chemotherapeutics cannot use without specialist doctor's advice.

The time and dosage given by the specialist doctor cannot be used. It cannot be used with other devices unless recommended by a specialist.

The oxygen concentrator is not used as a life support unit.

2. SAFETY GUIDELINES

- The device must only be used for purposes mentioned in this manual!
- Do not use this device as Life Supporting or Life Sustaining unit. Use your device according to your doctor's prescription and advices only, otherwise it may present a risk to the health of you (patient).
- Use of highly concentrated oxygen could be dangerous if not prescribed by authorized person. Please keep spare oxygen sources such as oxygen tube or bag ready for the patient considering possible malfunctioning of the device may occur anytime. Keep the device up right position during transportation to avoid damage in the cabinet. Switch off the device when not used by the patient.
- Working Conditions : Temperature 10 - 40 °C , Relative Humidity %20 - %60 RH , Altitude max 3000 m
- Storage Conditions : Temperature 10 - 60 °C , Relative Humidity %10 - %95 RH , Altitude max 4000 m

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- During transportation or when not in use unplug the power supply cord. Do not allow unauthorized person to open the cabinet of the device. If the electric supply line suffers from frequent voltage fluctuation, we suggest using Voltage Regulator for the device to avoid unwanted damages in the sensitive electronic parts of the device.
 - Do not unplug the power cable by hanging it on the device or by pulling the cable.
 - Do not leak water on the device or immerse the device into water. If you think there is a chance of water may come in contact with the device, please unplug the power cable and inform technical service.
 - Do not open the cabinet of the device when power cable connected. Only authorized technical service personnel is permitted to open the device. There are no parts inside the device which you can repair.
 - Do not put the nasal cannula under heavy material like bed or cushion.
 - Keep the device under surveillance when used by aged, child or disabled patient who may not able to inform any malfunctioning of the device may take place during the treatment.
 - Call your doctor immediately in case of slightest disturbance seen in the patient during treatment.
 - HEALTHTIME device produce highly concentrated oxygen which is not flammable but accelerates the combustion of material. To avoid all risk of fire do not smoke nearby the device and other accessories. HEALTHTIME should be kept away from all flames, incandescent sources and sources of heat, as well as combustible products like oil, grease, solvents, aerosols, etc.
 - The manufacturer or seller of the device is not responsible for the problems that may occur due to wrong, unnecessary or unintended use of the device or use of wrong electrical voltage which not comply with electrical safety rules.
 - Calibration requirement in the device during use of the device: During use of the device every 2000 hours or once in a year following calibration controls to be done by manufacturer. Oxygen Concentrator percentage b. Oxygen Flow rate c. Oxygen outlet pressure
 - Flow Cut-off Valve Against Fire: This feature available as standard in our device has double protection facility. The thermostat located in the compressor detect the increase in Motor temperature and shutdown the motor when the temperature reach to a deadline and a second thermostat located inside the device continuously measure the temperature inside the device and if any fire risk detected immediately stop the device and gives an alarm.
 - Children, Lactating Women and Pregnant must use the oxygen concentrator under the control of specialist doctor. Oxygen Concentrator use requires a physician's prescription and are not intended for life support use. Although oxygen therapy can be prescribed for patients of all ages, the typical oxygen therapy patient is older than 65 years of age and suffers from chronic obstructive pulmonary disease (COPD). Patients typically have good cognitive abilities and must be able to communicate easily. If the user is unable to communicate discomfort or unable to read and understand the concentrator labeling and instructions for use, then use is recommended only under the supervision of one who can. If any discomfort is felt while using the concentrator, patients are advised to contact their doctors. Patients are also advised to have back-up oxygen available in the event of a power outage or concentrator failure. .Device lifetime is 5 years. .Should be disposed according to weee (Waste Electrical and Electronic Equipment) Directive
- Brand of water-based lotions or ointments that are compatible with oxygen before and during oxygen therapy. Do not use petroleum or oil lotions or ointments for fire or burn vaccination. To avoid the risk of fire and burns, do not lubricate the oxygen concentrator fittings, fittings, pipes or other accessories.
- Use of this device at a relative humidity higher than the maximum declaration height or outside the temperature range of the statement temperature or at a relative humidity higher than the maximum statement relative humidity; it is expected to affect the flow rate and oxygen percentage and consequently adversely affect the quality of treatment.
- From children against choking risks caused by the cable and hose of the device, keep away from unconscious patients
- The device contains small and swallowable parts. Against the risk of swallowing keep away from children and unconscious patients
- If the device and its additional parts are used outside of their use purpose the manufacturer is not responsible for the faults that may occur
- After the use of spare parts other than the parts specified by the manufacturer the manufacturer is not responsible for the faults that may occur
- Disassemble or replacement of product parts not described in the user manual, from product safety damaged in case of combining with different products and the manufacturer is not responsible for any faults that may occur
- Keep the device away from children and pets
- The device is not suitable for use by more than one patient using the same parts.
- It should be used with new single-use components in different patients

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
3. INSTALLATION / USE

Please control the packing and the device on arrival. If there is any deformation or damage in the packing or in the inside materials, please inform the courier or the seller of the goods immediately. Take out from the package all components carefully controlling crash, breakage or any other damages.

NOTE: If Oxygen concentrator will not be used immediately, we suggest to keeping the components inside the package until use.

Let stand four (4) hours in user environment before start using the device.

Fill the water bowel with boiled and cooled water and close the lid. Take care of not over filling the water bowel out of its maximum mark. The water bowel must be cleaned with soap & water and rinsed properly once in a week. The water label of the bowel should control frequently and do not permit to go down to its minimum mark by adding water regularly.

Connect the intermediate cannula between device and patient unit (flowmeter unit). Place the water bowel in its position in patient' s unit as shown in picture. Plug in the electrical cable. Start the device by pushing the OFF-ON switch to ON position. Adjust the flow rate of concentrate oxygen with the help of inbuilt flow meter within the range of 0-5 lt./min. as per doctor's prescription. The ball located inside the flow meter will indicate the adjusted flow rate of oxygen. There is no matter that requires your intervention on the LCD screen of the device.

The total operating time and malfunctions of the device are directly displayed on this LCD screen. Please consult with "Troubleshooting" section of this manual if the device is not running properly.

To use a new device for the first time it is recommended to run the device for 30 min.

without patient. This period is necessary to achieve 90% oxygen concentration level at the initiation. Repeat the same procedure if the device is not used for a time period of 15 days or more. For daily use there is no need of following this instruction.

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Please the device at the most suitable place in the room or near to the patient. The device can be moved easily with the help of wheel fitted with the cabinet. Do not cover the concentrator with any covering material. There should be at least 30 cm. distances between the device & the furniture in the room and the walls to avoid blocking in environmental air inlet to the device. Do not block the air inlet to oxygen concentrator by any means. Do not place the device on soft floors like bed, sofa etc. it

may block air flows to the device. Keep the device always vertical position and do not use the device during bath.

Do not smoke near the device when running. Keep the device away from sources of fire and heat. It is dangerous to use the device with dirty filter or without filter. Never run the device with wet filter. Control the nasal cannula in a regular interval against twisting, leaking or bursting. Replace the damaged cannula with new one.

Like most other electronic devices too is very sensitive towards radio frequency waves (parasites), please be aware of using mobile phones when running the device. To get best efficiency from the device we suggest not shutting down and opening the device frequently in a short interval which may cause loss in life period of

the device too.

4. PREVENTIVE MAINTENANCE

Please switch off the device and unplug the power supply cable before maintenance operation start.

Do not use any liquid cleaning agent directly on the device. Petroleum based solvent & cleaner not to be used. Wipe the outer face of the cabinet with wet sponge and dry properly once in 30 days.

CLEANING OF DUST FILTER

Take out the dust filter placed at rear side of the device once in every 7 days.

The black color dust filter may clean from dust by simply shaking it and then wash with water only. Put in place after drying properly.

Replace the filter when it is torn or every 90 days.


CHANGING OF PALL FILTER

Replace the Pall filter with a new one every 60 days.

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

CHANGING OF NASAL CANNULA

Replace the Nasal Cannula with a new one every 30 days.

CLEANING AND CHANGING OF HUMUDIFIER

In every 7 days clean with soap and water , rinse properly and also Replace with a new one every 90 days

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
5. TROUBLE SHOOTING

OBSERVATION	ALARM SOUND	SCREENSHOT	PROBABLE REASON/CAUSES	SOLUTION
On opening there is no alarm and device is not running.	OK	NO IMAGE	Electrical power supply cord is not properly connected either with device or in plug. There is no electricity in the plug, fuse burnt out.	Connect the electrical cord properly both ends. Test if there is electric in the plug or put the cord in another plug. Change the fuse.
Alarm due to high temperature.	OK	HIGH TEMPERATURE	System may be turned to high temperature protecting mode.	Clean the black filter located at the rear side of the device. Check whether the ventilation is sufficient all around the device. Cooling fan may not working. Call authorized technical service. Keeps the device waiting for 40-50 mins. in switched off position then reopen, if problem continue call technical service.
High pressure alarm.	OK	HIGH PRESSURE	System may turn to high pressure protecting mode.	Control the Nasal Cannula against possible breaking, twisting. Control the water bowel against clogging. Check whether flow meter is open or closed. If the situations still continue please call Authorized Technical Service.
Low Pressure Alarm.	OK	LOW PRESSURE	System may turn to low pressure protecting mode.	Control the dust filter, if require clean it. Flow meter may set at 5 lt./min. reset it to a lower value. Please check, intermediate connecting hose may not be connected. If the problems still continue call technical service.
Low Oxygen Alarm	OK	LOW OXYGEN	A low oxygen condition may have occurred in the system.	Control the dust filter, pall filter, nazal cannula and flowmeter. If the problems still continue call technical service.
Call technical service if the suggested troubleshooters are not helpful for solution of problems. For problems not mentioned in this table, please consult with authorized technical service.				

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DO NOT STORE OR OPERATE THE DEVICE OUTSIDE THE SPECIFIED TEMPERATURE AND HUMIDITY CONDITIONS.

8. ACCESSORIES – PACKAGE INCLUDED

Package Included:

- 1 pc Oxygen Concentrator
- 1 pc Oxygen Cannula
- 1 pc Humudifier
- 1 pc Power Cord

Accessories: Cannula and Water container (humidifier)

Contact the authorized technical services for procurement.

Unsuitable parts or accessories may cause poor performance or malfunction of the device.

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9. ELECTROMAGNETIC COMPATIBILITY (EMC) DECLARATION

This equipment generates, uses and can radiate radio frequency (RF) energy. This equipment may cause electromagnetic interference if it is not installed and used in accordance with the manual. This device has been tested for Medical Products in accordance with EN-60601-1-2:2015 and has been determined to comply with acceptable limits. These limits indicate that the device provides an acceptable level of protection against electromagnetic interference (EMC) if the device is used as specified in the manual. This device may be affected by portable and mobile RF communication devices. This equipment must not be stored together with other equipment. See the following tables for more information about this device and EMC.

GUIDANCE AND MANUFACTURER'S DECLARATION – ELECTROMAGNETIC EMISSIONS

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment

Emisyon testi	Uyumluluk	Elektromanyetik ortam - kılavuz
Emissions RF CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. The device is suitable for use in all establishments including domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Emissions RF CISPR 11	Class B	
Harmonic emissions used for domestic purposes. IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

GUIDANCE AND MANUFACTURER'S DECLARATION – ELECTROMAGNETIC EMISSIONS

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment

Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharges IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30
Electrical fast transient/burst IEC 61000-4-4	Power supply lines: ± 2 kV Input/output lines: ± 1 kV 100 kHz repetition frequency	Power supply lines: ± 2 kV Input/output lines: ± 1 kV 100 kHz repetition frequency	Mains power quality should be that of a typical commercial or hospital environment.
Surges Line-to-line IEC 61000-4-5	± 0,5 kV, ± 1 kV	± 0,5 kV, ± 1 kV	Mains power quality should be that of a typical commercial or hospital environment.
Surges Line-to-ground IEC 61000-4-5	± 0,5 kV, ± 1 kV, ± 2 kV	± 0,5 kV, ± 1 kV, ± 2 kV	
Voltage dips IEC 61000-4-11	0% UT; 0.5 cycle 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°	0% UT; 0.5 cycle 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°	Mains power quality should be that of a typical commercial or hospital environment.
	0% UT; 1 cycle and 70% UT; 25/30 cycle Single phase: 0°	0% UT; 1 cycle and 70% UT; 25/30 cycle Single phase: 0°	
Voltage interruptions IEC 61000-4-11	0% UT; 250/300 cycle	0% UT; 250/300 cycle	
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Not: UT is the a.c. mains voltage prior to application of the test level.

Conducted RF IEC 61000-4-6	150kHz to 80MHz: 3Vrms 6 Vrms (in ISM and amateur radio bands) 80% Am at 1kHz	150kHz to 80MHz: 3Vrms 6 Vrms (in ISM and amateur radio bands) 80% Am at 1kHz	Portable and mobile RF communications equipments should be used no closer to any part of the device including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 0.35\sqrt{P}$ $d = 1.2\sqrt{80 \text{ MHz ile } 800 \text{ MHz}}$ $d = 1.2\sqrt{800 \text{ MHz ile } 2.7 \text{ GHz}}$ $d = 2.3\sqrt{P}$ Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	10V/m 80% Am at 1kHz	10V/m 80% Am at 1kHz	

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Note 1: 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.

b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

c. Calibration for current injection clamps shall be performed in a 150 Ω system.

d. If the frequency stepping skips over an ISM or amateur band, as applicable, an additional test frequency shall be used in the ISM or amateur radio band. This applies to each ISM and amateur radio band within the specified frequency range.

e. Applicable to ME EQUIPMENT and ME SYSTEMS with RATED input current less than or equal to 16 A / phase and ME EQUIPMENT and ME SYSTEMS with RATED input current greater than 16 A / phase.

Recommended separation distances between portable and mobile RF communications equipment and the device

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter(m)			
	150 kHz to 80 MHz (out ISM and amateur radio bands) $d = 1.2 \sqrt{P}$	150 kHz to 80 MHz (out ISM and amateur radio bands) $d = 0.6 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.7 GHz $d = 2.3 \sqrt{P}$
0,01	0,12	0,06	0,12	0,23
0,1	0,38	0,19	0,38	0,73
1	1,2	0,6	1,2	2,3
10	3,8	1,9	3,8	7,3
100	12	6	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: The ISM (industrial, scientific and medical) bands between 0,15 MHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.

Note 3: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Radiated RF	Test Frequency (MHz)	Band a) (MHz)	Service a)	Modulation b)	Modulation b) (W)	Distance (m)	Immunity Test Level (V/m)
IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipments)	385	380-390	TETRA 400	Pulse modulation b) 18 Hz	1,8	0,3	27
	450	380-390	GMRS 460, FRS 460	FM c) ± 5 kHz deviation 1 kHz sine	2	0,3	28
	710	704-787	LTE BAND1 13, 17	Pulse modulation b) 217 Hz	0,2	0,3	9
	745						
	780						
	810	800-960	GSM 800/900, TETRA 800, iDEN 820, LTE Bant 5	Pulse modulation b) 18 Hz	2	0,3	28
	870						
	930						
	1720	1700-1990	GSM 1800, CDMA 1900, GSM 1900, DECT, LTE Bant 1, 3, 4, 25, UMTS	Pulse modulation b) 217 Hz	2	0,3	28
	1845						
	1970						
	2450	2400-2570	Bluetooth, WLAN 802.11b/g/n, RFID 2450, LTE Bant 7	Pulse modulation b) 217 Hz	2	0,3	28
	5240	5100-5800	WLAN 802.11a/n	Pulse modulation b) 217 Hz	0,2	0,3	9
	5240						
	5785						

Note: If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.

a. For some services, only the uplink frequencies are included.

b. The carrier shall be modulated using a 50 % duty cycle square wave signal.

c. As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

The MANUFACTURER should consider reducing the minimum separation distance, based on RISK MANAGEMENT, and using higher IMMUNITY TEST LEVELS that are appropriate for the reduced minimum separation distance. Minimum separation distances for higher IMMUNITY TEST LEVELS shall be calculated using the following equation:

$$E = 6 / d \sqrt{P}$$

Where P is the maximum power in W, d is the minimum separation distance in m, and E is the IMMUNITY TEST LEVEL in V/m.

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www.elmaslarmedikal.com.tr

10. GUARANTEE (TERMS & CONDITIONS)

- 1- Guarantee period commences from the date of deliver of the product to client.
- 2- In case of malfunctioning within the warranty period, repair time is added to the warranty period. Maximum repairing period is 30 (Thirty) working days.
- 3- If you have any objection regarding the device, please return the device to the company from where you bought the device or to the manufacturer by prepaid cargo along with the sales certificate and the guarantee certificate.
- 4- Problems appeared due to wrong or misuse of the device is not covered under warranty.
- 5- Opening of the device by unauthorized person will cause device to remain outside of the scope of the warranty.
- 6- The guarantee for the device will only be valid once the consumer accepts after reading all terms and conditions mentioned here.
8. Consumable parts are not covered under warranty.

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12. MANUFACTURER - AUTHORIZED TECHNICAL SERVICES

ELMASLAR A.Ş.

ASO 1.OSB BABÜRŞAH CD NO : 17 SİNCAN - ANKARA / TURKEY

Tel: +90 312 394 80 01 Faks: +90 312 394 80 04

info@elmaslarmedikal.com.tr

www.elmaslarmedikal.com.tr

ELMASLAR A.S.

ASO 1.OSB BABURSAH CD. NO: 17 SINCAN / ANKARA / TURKEY

Tel: +90 312 394 80 01 Fax: +90 312 394 80 04

info@elmaslarmedikal.com.tr • elmaslarmedical@gmail.com

www.elmaslarmedikal.com.tr