


555



888



FLEX 5



TEKNİK ÖZELLİKLER

HORTUM TİPİ	555	888	FLEX 5
STERİLİZASYON SICAKLIĞI	< 70°C	< 134°C (OTOKLAVLANABİLİR)	< 70°C
HAMMADDE	HORTUM: POLYETHYLENE, BAŞLIK: ELASTOMER	HORTUM: COPOLYETER, BAŞLIK: SİLİKON	HORTUM: POLYETHYLENE, BAŞLIK: ELASTOMER
ÇAP	HORTUM: 19 mm, BAŞLIK: 22 mm	HORTUM: 19 mm, BAŞLIK: 22 mm	HORTUM: 19 mm, BAŞLIK: 22 mm
UZUNLUK	STANDART 180 cm 0 - 300 mm OPSİYONEL	STANDART 180 cm 0 - 300 mm OPSİYONEL	STANDART 180 cm 0 - 300 mm OPSİYONEL
KULLANIM ALANLARI	UYKU APNESİ UYGULAMALARI	ANESTEZİ ULTRASONİK NEBULİZATÖR UYKU APNESİ UYGULAMALARI	UYKU APNESİ UYGULAMALARI