



# SLEEPTIME+ BPAP/BPAP ST/BPAP ST-30



Product Page

The purpose in using this device is to keep the upper airways open during sleep by applying continuous and constant air pressure.



### Wide LCD Screen

Thanks to its wide LCD screen, all parameters can be monitored from a single screen.



### SD Card Assisted

All data is stored on the SD card and easy reporting is provided for the desired periods.



### Free User Interface Design

The user interface is available free of charge from the site.



### Auto Start/Stop Function

Automatic Start/Stop can be made thanks to the breath sensitive sensor.



### RESlex Pressure Adjustment

When the patient's constant pressure changes, the device automatically detects and delivers the correct pressures, increasing patient compliance.

### Technical Specifications

Product Model	BPAP S	BPAP-ST
Working Mode	CPAP/BPAP S	CPAP/BPAP S/ BPAP-ST
Automatic Pressure		No
RESlex		Yes
Altitude Compensation		Yes
Ti Control		Yes
Working Pressure (hPa) Range	4-25	4-25 / 4-30
Max. Flow		≤ 40
SD Card Capacity		4 gb
Noise Level		< 30dB
Max. Delivery Temperature		< 43°C
Dimensions		140×190×115mm 280×190×130mm
Weight		2.5 kg (with Humidifier)
AC Power Consumption		100 – 240 V AC 50/60 Hz, 3.0 A Max
Protection against Electric Shock		CLASS II
Standard Compliance		IEC 60601-1 ISO 17510
Stock Code	N1221004	N1221006