



NEBTIME CN02



Product Page

This compact device was created to effectively deliver the medication prescribed by the doctor to the patient.



Compact Design

Thanks to its ergonomic design, it provides easy of use and storage.



Converting the Prescribed Medication

Thanks to its powerful motor, it delivers the prescribed medication to the patient within minutes.



Oil Free Piston Pump

Its special compressor does not require maintenance and allows use for many years.



Optimum Time Option

It ensures that the prescribed dose is given to the patient at the right time.



Handle for Easy Use

Thanks to its handle, it is very easy to carry.



Storage for Accessories

Thanks to its storage, all consumables are kept together hygienically.

Technical Specifications

| | |
|-----------------------------|--------------------------|
| Voltage | 230 VAC \pm 10%, 50 Hz |
| Power | 100W |
| Medication Capacity | 5 ml |
| Particle Size | 0,5 – 10 μ m |
| MMAD | 4 μ m |
| Noise Level Below | < 60 dB |
| Maximum Medication Capacity | 13 ml (cc) |
| Average Nebulization Rate | min 0,2 ml / min |
| Compressor Pressure Range | 30 ~ 36 psi (max)/bar |
| Operation Pressure Range | 8 ~ 16 psi/bar |
| Liter Flow Range | 6 ~ 8 lpm |
| Operating Temperature Range | 10°C ~ 40°C |
| Operating Humidity Range | 10 ~ 90 % RH |
| Device Dimensions | 32x17x11 cm |
| Box Dimensions | 34x19x12 cm |
| Stock Code | NB02B000 |
| Warranty Period | 2 Years |
| Power Consumption | 0,70 A |
| Weight of Device | 1.8 kg |
| Gross Weight | 2 kg |