



Optional stand

NEBTIME CN03



Product Page

This compact device was created to effectively deliver the medication prescribed by the doctor to the patient.

- 24 Hours Continuous Use**
The product was created with hospitals in mind, and it may be used in a medical environment for up to 24 hours in a row.
- Outlet Pressure up to 6 Bar**
It provides pressure up to 6 bar thanks to its adjusting knob for patients who need high pressure.
- Medicine Transfer to 3 Patients**
If needed, 1 device with its specially designed apparatus can be used on 3 patients without pressure loss.
- Low Cost**
Since 1 device allows the use of 3 patients, the cost of the device per patient is low.

Technical Specifications

Compressor Type	Oilless Piston Pump
Flow Range	30 lt/min \pm %10
Pressure Range	3,0 bar \pm %10
Nebulization Range	0,50 ml/min \pm %10
Noise Level	60 dB \pm %10
Particle Size	0,5 - 10 μ m
MMAD	3 μ m \pm %10
Voltage	230 VAC \pm %10, 50 Hz/60Hz
Power	100 W
Operation Type	Continuous
Operating Temperature Range	10°C ~ 40°C
Device Dimensions	29x17x25 cm
Box Dimensions	35x20x27 cm
Gross Weight	5 kg
Stock Code	NB03B000
Warranty Period	2 Years