



NEBTIME CN04



Product Page

This compact device was created to effectively deliver the medication prescribed by the doctor to the patient.



Compact Design

Thanks to its ergonomic design, it provides easy of use and storage.



Converting the Prescribed Medication

Thanks to its powerful motor, it delivers the prescribed medication to the patient within minutes.



Oil Free Piston Pump

Its special compressor does not require maintenance and allows use for many years.



Optimum Time Option

It ensures that the prescribed dose is given to the patient at the right time.



Handle for Easy Use

Thanks to its handle, it is very easy to carry.



Storage for Accessories

Thanks to its storage, all consumables are kept together hygienically.

Technical Specifications

Voltage	230 VAC \pm %10, 50 Hz
Power	100W
Medication Capacity	5 ml
Particle Size	0,5 – 10 μ m
MMAD	4 μ m
Noise Level Below	< 60 dB
Maximum Medication Capacity	13 ml (cc)
Average Nebulization Rate	min 0,2 ml / min
Compressor Pressure Range	30 ~ 36 psi (max.)/bar
Operation Pressure Range	8 ~ 16 psi/bar
Liter Flow Range	6 ~ 8 lpm
Operating Temperature Range	10°C ~ 40°C
Operating Humidity Range	10 ~ 90 % RH
Device Dimensions	29x16x10 cm
Box Dimensions	32x17x11 cm
Stock Code	NB04B000
Warranty Period	2 Years
Power Consumption	0,70 A
Weight of Device	1.5 kg
Gross Weight	1.7 kg