



SLEEPTIME+ CPAP/AUTO CPAP



Product Page

The purpose in using this device is to keep the upper airways open during sleep by applying continuous and constant air pressure.



Wide LCD Screen

Thanks to its wide LCD screen, all parameters can be monitored from a single screen.



SD Card Assisted

All data is stored on the SD card and easy reporting is provided for the desired periods.



Free User Interface Design

The user interface is available free of charge from the site.



Auto Start/Stop Function

Automatic Start/Stop can be made thanks to the breath sensitive sensor.



RESlex Pressure Adjustment

When the patient's constant pressure changes, the device automatically detects and delivers the correct pressures, increasing patient compliance.

Technical Specifications

Product Model	CPAP	AUTO CPAP
Working Mode	CPAP	CPAP/AUTO
Automatic Pressure	No	Yes
RESlex		Yes
Altitude Compensation		Yes
Working Pressure (hPa) Range		4-20
SD Card Capacity		4 gb
Noise Level		< 30dB
Max. Delivery Temperature		< 43°C
Dimensions		140×190×115mm 280×190×130mm
Weight		2.5 kg (with humidifier)
AC Power Consumption		100 – 240 V AC 50/60 Hz, 3.0 A Max
Protection against Electric Shock		CLASS II
Standard Compliance		IEC 60601-1 ISO 17510
Stock Code	N1221001	N1221003