



HOSES



Product Page

It provides superior performance and maximum comfort. Lighter design and effective treatment it offers you and your family a quality night's sleep.

Technical Specifications

Hose Type	555	888	FLEX 5
Sterilization Temperature	< 70°C	< 134°C (Autoclavable)	< 70°C
Raw Material	Hose: Polyethylene Nozzle: Elastomer	Hose: Copolyester Nozzle: Silikon	Hose: Polyethylene Nozzle: Elastomer
Diameter	Hose: 19mm Nozzle: 22mm	Hose: 19mm Nozzle: 22mm	Hose: 19mm Nozzle: 22mm
Length	Standard 180 cm 0-300 mm Optional	Standard 180 cm 0-300 mm Optional	Standard 180 cm 0-300 mm Optional
Fields of Usage	Sleep Apnea Treatment	Anesthesia Ultrasonic Nebulizer Sleep Apnea Treatment	Sleep Apnea Treatment