



# HEALTHTIME OC10



Product Page

The Oxygen Concentrator is a device that produces oxygen with a purity of  $93 \pm 3\%$  by using environment air containing 21% oxygen.

**15M**

### Unique Patient Unit Design

It increases patient comfort without any loss of pressure and purity, with its unit extended up to 15 meters. (Optional)



### Ergonomic Design

It provides ergonomic use with its hidden handle design and wheels made of soft material.



### Lightness

It is one step ahead of other oxygen concentrators with a weight of only 19 KG.



### Noise Level

It has the lowest noise level on the market, with the contribution of its unique extendable patient unit.



### Polycarbonate Outer Case

Unlike other manufacturers, we use polycarbonate to remove harmful gases and bad odor.

### Technical Specifications

Oxygen Flow Range	1~ 10l / min
Oxygen Concentration	%93 ±3
Maximum Outlet Pressure	9 PSI ±1
Voltage	230 VAC ±%10, 50 Hz
Current	2,60 AMP ±%10 500W
Operating Temperature Range	10°C ~ 40°C
Operating Humidity Range	10 ~ 90 % RH
Weight Net/Gross	19 kg/22 KG
Noise Level	< 45 Db
Device Dimensions	28x34x52 cm
Box Dimensions	38x50x61 cm
Stock Code	OC01B003
Warranty Period	2 Years
Digital Display	Yes
Power Consumption	500W