



HEALTHTIME OC5



Product Page

The Oxygen Concentrator is a device that produces oxygen with a purity of $93\pm 3\%$ by using environment air containing 21% oxygen.

15M

Unique Patient Unit Design

It increases patient comfort without any loss of pressure and purity, with its unit extended up to 15 meters. (Optional)



Ergonomic Design

It provides ergonomic use with its hidden handle design and wheels made of soft material.



Lightness

It is one step ahead of other oxygen concentrators with a weight of only 14 KG.



Noise Level

It has the lowest noise level on the market, with the contribution of its unique extendable patient unit.



Polycarbonate Outer Case

Unlike other manufacturers, we use polycarbonate to remove harmful gases and bad odor.

Technical Specifications

Oxygen Flow Range	1~ 5l / min
Oxygen Concentration	%93 ±3
Maximum Outlet Pressure	8 PSI ±1
Voltage	230 VAC ±%10, 50 Hz
Current	1,40 AMP ±%10
Operating Temperature Range	10°C ~ 40°C
Operating Humidity Range	10 ~ 90 % RH
Weight Net/Gross	14 KG/18 kg
Noise Level	< 45 Db
Device Dimensions	28x34x52 cm
Box Dimensions	38x50x61 cm
Stock Code	OC01B002
Warranty Period	2 Years
Digital Display	Yes
Power Consumption	370W