



5L

OXYTIME Q5



Product Page

The Oxygen Concentrator is a device that produces oxygen with a purity of $93\pm 3\%$ by using environment air containing 21% oxygen.



Digital Display

Thanks to the digital display, it provides quick access to all information and facilitates control with alarm and timer.



Ergonomic Design

Designed with an easy-to-use handle, Oxytime Q5 offers an ergonomic experience.



Lightness

It is one step ahead of other oxygen concentrators with a weight of only 14 KG.



Lockable Castor Wheels

It's not only portable, but it's also exceptionally safe to use thanks to its four lockable castor wheels.

Technical Specifications

Oxygen Flow Range	1~ 5l / min
Oxygen Concentration	%93 \pm 3
Maximum Outlet Pressure	8 PSI \pm 1
Voltage	230 VAC \pm 10, 50 Hz
Current	1,60 AMP \pm 10
Operating Temperature Range	10°C ~ 40°C
Operating Humidity Range	10 ~ 90 % RH
Weight Net/Gross	14 KG/16 kg
Noise Level	< 45 Db
Device Dimensions	37x21x58 cm
Box Dimensions	45x36x65 cm
Stock Code	OC01B004
Warranty Period	2 Years
Digital Display	Yes
Power Consumption	370W