



20-30  
LPM

## LIFETIME SA01PB/SA02PB



Product Page

Ideal for ambulances, emergency services, examination rooms, inpatient services, outpatient clinics and home care.



### Analog Display

Its analog screen provides the most accurate way to monitor the adjusted pressure with a sensitive vacuum meter.



### Adjustment Knob

You can reach the desired pressure in the shortest time with the sensitive adjustment knob.



### Li-ion Battery

Provides up to 60 minutes of use on a single full charge. No need for battery replacement for a long time.



### Wall Bracket

Thanks to its wall bracket that comes as a standard in clinics and ambulances, stable use is ensured.



### LED Charge Indicator

Thanks to the LED charge indicator, the battery charge rate can be monitored continuously in 5 stages.



### Accessory Box

The accessory box provides easy transportation and use of all consumables with the device, especially for mobile uses.

### Technical Specifications

Product Model	SA01PB	SA02PB
Flow Rate	20 LPM $\pm$ %10	30 LPM $\pm$ %10
Noise Level	<50 dB	<60 dB
Weight	4 kg	4.5 kg
Battery	90 min	60 min
Max. Vacuum	-650 $\pm$ 5 mmHg	
Power Consumption	30W	60W
Compressor Type	Oilless Piston Pump	
Battery Type	Li-ion	
Voltage	12 VDC	
Jar Capacity	1 lt	
Box Dimensions	40x20x33 cm	
Gross Weight	5 kg	6 kg
Warranty Period	2 years	
Stock Code	PB01B000	PB02B000