



LIFETIME SA01PT



Product Page

Is a suction unit used for body fluids, oral, nasal and tracheal aspirations in adults and children.



Compact Design

Thanks to its ergonomic design, it provides ease of use and storage.



Adjustment Knob

You can reach the desired pressure in the shortest time with the sensitive adjusting knob.



Autoclavable Jar

Thanks to its autoclavable polycarbonate jars at 121 degrees, it provides hygienic use for many times.



Oil Free Piston Pump

Its special compressor does not require maintenance and allows use for many years.

Technical Specifications

Compressor	Oilless Piston Pump
Voltage	230 VAC \pm %10 , 50 Hz
Max. Free Air Flow	12 lt / min \pm %10
Max. Vacuum	-550 mmHg \pm %5
Power Consumption	0,60 A
Vacuum Setting	Yes
Collection Jar	1 lt
Float System	Yes
Hydrophobic Filter	Yes
Standard	ISO 10079 –1/ CE
Device Dimensions	20x15x25 cm
Weight of Device	2.5 kg
Box Dimensions	34x16x28 cm
Gross Weight (Device)	2.8 kg
Stock Code	PT01B000
Warranty Period	2 years