



LIFETIME SA02PT/SA03PT



Product Page

Designed to evacuate the fluid accumulated in the operated areas of the patients during the surgical operations.



Easy Operation

It is ready for use in a short time with one-button opening and sensitive adjustment knob.



Suitable for Clinical Use

Apart from home use, it is also suitable for use in clinics thanks to its high vacuum power.



Adjustment Knob

You can reach the desired pressure in the shortest time with the sensitive adjustment knob.



Ergonomic Design

Thanks to its compact design and handle, it provides ease of transport at home and in clinics.

Technical Specifications

Product Model	SA02PT	SA03PT
Compressor	Oilless Piston Pump	
Voltage	230 VAC \pm %10 , 50 Hz/60Hz	
Max. Vacuum	-640 mmHg \pm %5	
Max. Free Air Flow	15 lt/min \pm %10	30 lt/min \pm %10
Vacuum Setting	Yes	
Collection Jar	1-2 lt	1-2 lt
Vacuum Entrance	Single	Single
Float System	Yes	
Hydrophobic Filter	Yes	
Standard	ISO 10079 – 1 / CE	
Device Dimensions	29x17x25 cm	
Weight of Device	4.2 kg	6.3 kg
Warranty Period	2 years	
Stock Code	PT02B000	PT03B000