



Optional stand

NEBTIME UN300A



Product Page

Evaporation of liquids and medicine in solution form allows it to be used in the treatment of patients.



Lightness

It is one step ahead of other ultrasonic nebulizers with a weight of only 1.2 KG.



Noise Level

The Nebtime UN300A is the device with the lowest noise level on the market.



Home/Optional Hospital Use

Designed with hospital staff in mind, Nebtime UN300A offers a comfortable use for the hospital as well.



Smooth Design

Elmaslar Medical Devices has won numerous design awards for its ergonomic design approach, which is evident in all of its products.

Technical Specifications

Voltage	230 VAC±%10, 50 Hz/60Hz
Power	50W
Crystal Operating Frequency	1,68 MHz ±%5
Particle Size	0,5 – 6 µm
Average Nebulization Rate	0-5 ml/dk ±%10
Aerosol Tank	300 ml
Weight	1.2 kg
Gross Weight	1.5 kg
Operating Temperature Range	10°C ~ 40°C
Storage Temperature Range	10°C ~ 40°C
Power Consumption	50W
Operation Type	Continuous
Device Dimensions	26x16x20 cm
Box Dimensions	36x18x23 cm
Stock Code	UN02B000
Warranty Period	2 Years
Noise Level	< 40 dB
Medicine Cup Capacity	70 ml
Maximum Medicine Capacity	30 ml